INSTRUCTIONS FOR USE (ENGLISH)
MANUAL DE INSTRUCCIONES DE USO (ESPAÑOL)

ESCAPE
Aspiration Catheter

STERILE EO

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
- Do not resterilise
- No reesterilizar

- Do not use if packaging is damaged
- No utilizar si el embase está dañado

- Keep away from sunlight
- Mantener lejos del sol

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CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

STERILE. Sterilised with ethylene oxide gas. Non pyrogenic. For one procedure only. Do not re-sterilise. Do not use opened or damaged packages. Destroy product after use. Store in a dry place below 40° C, keep away from light. Refer to accompanying Instructions for Use.

DESCRIPTION

The Escape Aspiration Catheter is a single-user design, dual lumen catheter. It has a distal radiopaque tip marker and proximal luer-lock port, the proximal luer-lock port is for connection of the Aspiration Line (supplied) and Aspiration Syringe (supplied).

INDICATIONS AND INTENDED USE

The Escape Catheter is indicated for use in the central and peripheral circulatory system, including saphenous vein grafts to:

- Contain and aspirate embolic material (thrombus/debris) while performing percutaneous transluminal coronary angioplasty, percutaneous transluminal angioplasty and/or stenting procedures.

WARNINGS

- Prior to use, the packaging and product should be inspected for signs of damage. Never use damaged product or product from a damaged package.
- The Escape Catheter should be handled carefully. Prior to use inspect the Escape Catheter carefully for bends, kinks, or other damage. Do not use a damaged Escape Catheter.
- The Escape catheter is designed and intended for single use. Check that the sterile pouch is not damaged before use. Do not re-sterilise or re-use, destroy product after use. Re-use & re-sterilization may compromise the device performance.

PRECAUTIONS

- The Escape Aspiration Catheter and accessories should be used in conjunction with fluoroscopic guidance and proper anticoagulation agents.
- This device is designed and intended for single patient use only. DO NOT re-sterilise and/or re-use.
- As in any elective coronary intervention, it is recommended that the patient have a mean systolic blood pressure greater than or equal to 90 mm Hg in concomitant of IV pressors or Intra Aortic Balloon Pump augmentation.

As with most percutaneous interventions, other potential adverse events include: Death, Myocardial Infarction, Coronary or Bypass Graft Thrombosis or Occlusion, Myocardial Ischemia, Stroke/CVA, Emergent or Non-emergent Bypass Graft Surgery, Emboli (air, tissue or thrombotic), Dissection, Arterial Perforation, Arterial Rupture, Ventricular Fibrillation, Hemorrhage, Hypotension, Pseudoaneurysm at Access Site, Arteriovenous Fistula, Infection at Access Site, Other Bleeding Complications at Access Site.

MATERIALS REQUIRED FOR USE WITH THE ESCAPE ASPIRATION AND INFUSION SYSTEM

- Minimum 6F or 7F Arterial or Venous Introducer Sheath and minimum 6F or 7F Femoral Guiding Catheter in the appropriate configuration to cannulate the vessel (preferably with side holes if ostial narrowing is present or the guide catheter is occlusive).
- Push pull or Rotating Hemostatic Valve
- Heparinized normal saline.

Recommendation:

- If the use of a second guide wire in addition to the Escape Aspiration Catheter is required, please note that an increased size of guiding catheter should be used in order to avoid friction issues.
- Example: for a 6F Escape Aspiration Catheter, use a 7F Guiding Catheter.

DIRECTIONS FOR USE**Preparation of the Escape Catheter**

The Escape Aspiration Catheter is supplied with an Aspiration Line, one locking Aspiration Syringe and 2 filters.

1. Remove the Escape Catheter and accessories from the package.

2. Fill the Aspiration Syringe with approximately 5-10 ml of heparinized saline and attach the aspiration line and syringe to Escape Catheter.
3. Open the stopcock on the Aspiration Line and flush the entire length of the Escape Catheter using all of the heparinized saline contained in the Aspiration Syringe. Close the stopcock.
4. Verify that the stopcock on the Aspiration Line is in the closed position. Retract the plunger of the Aspiration Syringe and pull until it locks at the fully extended position. The Escape Catheter is completely prepped and is ready for use.
5. Perform aspiration using the Escape Catheter:
 - a. Load and advance the prepped Escape Catheter over the guide wire to the tip of the Guiding Catheter.
 - b. Under fluoroscopy advance the Escape Catheter and position the distal tip marker proximal above the embolic particles. Stop advancement of the Escape Catheter if any resistance is encountered.
 - c. Begin aspiration by opening the stopcock on the Aspiration Line. Slowly retract the Escape Catheter towards the Guiding Catheter. Blood will enter the Aspiration Syringe until the entire vacuum is gone (or the Aspiration Syringe is filled).

NOTES:

- Should blood not begin filling the syringe within 5 seconds, check the Guiding Catheter tip placement. Unseat the Guiding Catheter if necessary to resume flow.
 - If no blood is aspirated as a result of unseating the Guiding Catheter, turn the stopcock off and remove the Escape Catheter outside of the patient, either flush the aspiration lumen or use a new Escape Catheter. DO NOT flush the system while the catheter is still inside the patient vasculature.
 - During performance testing, the Escape Catheter demonstrated the ability to evacuate fluid and debris at a minimum rate of 1 ml/second.
- d. After completing the aspiration process turn the Aspiration Line stopcock to close off the Aspiration Syringe.
6. Remove the Escape Catheter
- a. Withdraw the Escape Catheter
- NOTE:* Especially for peripheral intervention (or non-coronary): the guide wire must be withdrawn simultaneously with the catheter to avoid any loop creation.
- b. Slowly retract and remove the Escape Catheter. If necessary, loosen the Tuohy Borst of the Haemostatic valve to allow easy withdrawal of the distal shaft.
- NOTE:* Remove the Aspiration Syringe and re-flush the aspiration lumen and wire lumen of the Escape Catheter with heparinized saline before each re-use Empty the Aspiration Syringe, reattach to extension line, close the stopcock, and retract the plunger to the fully extended lock position
7. Remove the catheters and follow standard hospital practice for management of the insertion site.

HOW SUPPLIED

The Escape Catheter is supplied with an Aspiration Line, one locking Aspiration Syringe and 2 filters; it is supplied sterile and non-pyrogenic in a tray in its own pouch. The Escape Aspiration Catheter including all of its components are intended to be used for one procedure only.

NOTE ALTHOUGH THE ESCAPE ASPIRATION CATHETER, HEREAFTER REFERRED TO AS 'PRODUCT', HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, DISA VASCULAR 2015 (PTY) LTD HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. DISA VASCULAR 2015 (PTY) LTD THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE. DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT. WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE NO PERSON HAS ANY AUTHORITY TO BIND DISA VASCULAR 2015 (PTY) LTD TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected

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The Escape Aspiration Catheter 6F / 7F is a product of DISA Vascular 2015 (Pty) Ltd
El Catéter de Aspiración ESCAPE 6F / 7F es un producto de DISA Vascular 2015 (Pty) Ltd