INSTRUCTIONS FOR USE
INSTRUCCIONES DE USO

CAPECROSSnc
R_x PTCA Balloon Catheter

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
- Keep away from sunlight
- Mantener lejos del sol

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


STERILE. Sterilised with Ethylene Oxide gas. Non-pyrogenic. Does not contain phthalates. For one procedure only. Do not use opened or damaged packages. Destroy product after use. Do not resterilise.

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

1. DESCRIPTION

The Cape Cross NC is a non-compliant Percutaneous Transluminal Coronary Angioplasty (PTCA) catheter. The distal shaft comprises two lumens. One is used for inflation of the balloon with contrast medium. The second lumen, in the distal shaft, permits the use of a guidewire (0.014" max.) to enable advancement of the catheter to and through the stenosis to be dilated. The balloon has two radiopaque markers to aid in positioning the balloon in the stenosis. The distal part of the Cape Cross NC is coated with hydrophilic coating, which is activated when wet, to facilitate access to the lesion site.

The proximal shaft consists of a stainless steel hypotube. Proximal visual markers located approximately 90cm and 100cm from the distal tip aid catheter positioning for brachial and femoral techniques.

2. INDICATIONS

The Cape Cross NC catheter is intended to be used to dilate diseased segment(s) in a coronary artery, or coronary bypass graft, to improve coronary luminal diameter and thereby improve myocardial perfusion. It is used in patients with clinical symptoms of myocardial ischaemia related to the pathological condition of one or more coronary arteries or grafts. It may be used to dilate a vessel alone (POBA), or may be used for balloon dilatation of a stent after implantation. More details are provided in the Instructions below (Section 6 and 7).

Note: Bench testing was conducted with the Cape Cross NC and the SolarFlex stent. Consideration should be taken when this device is used with different manufacturers' stents due to differences in stent design. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.

3. CONTRA-INDICATIONS

- Unprotected left main coronary artery,
- Coronary artery spasm in the absence of a significant stenosis

4. WARNINGS

- The Cape Cross NC catheter is designed and intended for single use. Carefully examine all equipment to be used during the procedure, including the catheter, to verify proper function. Verify that the catheter and sterile pouch have not been

damaged. Do not use if sterile package is damaged. Verify that the catheter size is suitable for the specific procedure for which it is intended.

- Do not resterilize or reuse; destroy product after use. Reuse & resterilization may compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.
- Use the product before the expiry date clearly shown on the packaging.
- PTCA procedures should only be performed in hospitals where emergency coronary artery bypass graft surgery can be immediately performed in the event of potentially injurious or life-threatening complications.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during the procedure, as treatment of this patient population carries special risk.
- When the PTCA balloon catheter is exposed to the vasculature, it should be manipulated while under high-quality fluoroscopic observation.
- If resistance is met during manipulation, determine the cause of the resistance before proceeding. Attempts to advance or retract the catheter under resistance may result in damage to the vessel and/or damage/separation of the catheter.
- Use only diluted contrast medium to inflate the balloon. Do not use air or any gaseous medium.
- To reduce the potential for damage in coronary vessels and to avoid any balloon rupture, the balloon inflation pressure should not exceed the rated burst pressure (RBP) (cf. compliance chart). The RBP is based on results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their RBP. The probability of balloon rupture increases if this pressure is exceeded. The compliance chart reflects diameters and pressures at or below which 99% of balloons of balloons will not burst (95% confidence).
- Inflation beyond RBP may also increase the resistance to withdraw the catheter from the guiding catheter due to flattening of the balloon when applying vacuum.
- Use of a pressure-monitoring device is recommended to prevent over-pressurization.
- Do not attempt to advance or retract the catheter unless the balloon is fully deflated under vacuum. Attempts to reposition a fully or partially inflated balloon may result in severe vessel damage.
- Do not use, or try to straighten, a bent or kinked catheter: not following these instructions could result in shaft rupture. In case of defective product, use another catheter.
- Treatment of moderately or heavily calcified lesions increases the risk of vessel closure, vessel trauma, balloon burst, balloon entrapment and associated complications.
- Treatment of patients with totally obstructed arteries, diffuse or multiple lesions or multi-vessel disease, thrombus inside the target vessel lumen or presenting in acute

cardiogenic shock should be carefully considered as these treatments may be associated with increased risk of complications.

- During withdrawal of the catheter, it is recommended that a saline-soaked gauze is held around the exposed catheter shaft while the catheter is pulled through the gauze to remove blood or any other residues. Prior to reinsertion, the balloon should be submerged in sterile, heparinised normal saline to reactivate the coating.
- Care should be taken when handling the distal part of the catheter to prevent damage and prematurely removing the balloon protector.
- Before withdrawing the balloon catheter, visually confirm complete balloon deflation by fluoroscopy. Wait at least 2 seconds after the contrast medium has exited the balloon to be sure the balloon has fully retracted.

5. PRECAUTIONS

- Note the “Use by” date specified on the package.
- These devices should be used only by physicians trained in PTCA.
- Ensure that the medical team is trained on PTCA to avoid any error in choosing equipment.
- Inspect the catheter prior to use for any kinks, curves or potential catheter damage, which could alter the catheter’s performances (see procedure for preparing the Cape Cross NC as specified below) and ensure that its size is suitable for the specific procedure for which it is to be used.
- Do not attempt to inflate the balloon before it is correctly positioned on the level of the lesion to treat.
- The Cape Cross NC catheter must be used in patients who have been provided with appropriate anticoagulant and coronary vasodilator therapy as needed.
- With 4.5 mm balloon dilatation catheters, some increased resistance may be noted upon insertion or withdrawal into or out of the guiding catheter. Choosing a larger guiding catheter size may minimize this.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn from the patient before removing any additional equipment.
- Cape Cross NC catheters 3.0mm or larger are not suitable for simultaneous use within a 6F guide catheter.
- Care should be taken not to over-tighten the haemostatic adaptor around the catheter shaft as lumen constriction may occur, affecting inflation and deflation of the balloon.
- Do not reinsert the Cape Cross NC catheter into the coil dispenser after procedural use.
- The design and construction of these catheters do not provide the user with distal pressure monitoring capability.

6. PREPARATION FOR USE

Prior to use, examine all equipment for defects. Do not use any defective equipment.

A. Catheter selection:

The inflation diameter of the balloon must not exceed the diameter of the coronary artery proximal and distal to the stenosis.

B. Catheter preparation:

Select an appropriate balloon catheter for the target vessel.

Remove the device from the sterile pouch.

Remove the catheter from the protective hoop. Use care when removing the catheter to avoid damage (e.g. shaft kink).

Remove the balloon protector and mandrel by grasping the catheter just proximal to the balloon catheter (at the proximal balloon catheter bond site). With the other hand, gently grasp the balloon protector and remove distally. Remove the mandrel distally after removing the balloon protector.

The guide wire lumen may be flushed by filling a 10 or 20cc syringe with normal sterile heparinised saline solution, attaching to it a flushing needle and inserting the needle into the distal tip of the catheter. Then flush the lumen by gently applying pressure to the syringe.

Prepare a 50% solution of contrast medium in sterile heparinised saline, in accordance with the contrast medium manufacturer's instructions.

Connect a stopcock to the hub of the Cape Cross NC catheter.

Attach a 10 or 20 cc luer lock syringe with 1 to 10 cc of a 50% solution of contrast medium in sterile heparinised saline, to the stopcock.

Position the syringe vertically with the plunger pointing upwards. Apply strong suction with the syringe to evacuate all air from the balloon.

Maintain the suction for 15 to 20 seconds and make sure that no more bubbles are seen passing through the diluted contrast medium.

Release the plunger carefully and, if necessary, repeat the procedure, in case of residual bubbles. Balloon purging is complete when the flow of bubbles stops. If the bubbles do not stop, inspect the system for possible leaks. In case of defects, do not use the catheter.

Maintain the vacuum in the Cape Cross NC catheter and close the stopcock. Disconnect the syringe.

Fill an inflation device with 3 to 10cc of a 50% solution of contrast medium in sterile heparinised saline according to the manufacturer's instructions.

Connect the inflation device, avoiding the presence of bubbles. Apply negative pressure and close the stopcock. Purge any air from the inflation device through the stopcock. Leave the catheter under negative pressure until it is ready to be used.

Note: The distal portion of the Cape Cross NC is covered with a hydrophilic coating. Immerse the distal part of the catheter in a sterile saline solution before use to activate the coating. Do not reinsert into dispenser.

CAUTION: All air must be removed from the balloon and displaced with contrast prior to inserting into the body, otherwise, complications may occur.

7. INSTRUCTIONS FOR USE

Attach the haemostasis valve (Y-connector) to the guiding catheter, which has been previously placed in the vasculature. The choice of the guidewire and the haemostasis valve remains the responsibility of the physician within a compatibility limit of 0.014" maximum for the guidewire.

Insert the guidewire into the guiding catheter. Under fluoroscopy, position the guidewire through the guiding catheter and across the stenosis according to accepted PTCA techniques.

To avoid any movement of the guidewire, tighten the knob of the haemostasis valve so that it closes firmly around the guidewire.

Insert the proximal end of the guidewire into the distal tip of the Cape Cross NC catheter. The guidewire will exit through the guidewire exit port located approximately 25cm proximal to the dilatation catheter tip.

Open the valve and advance the Cape Cross NC to the distal end of the guiding catheter. Caution: Care should be taken when inserting the Cape Cross NC catheter into the haemostasis valve in order to avoid kinking.

Two markers located on the proximal part of the shaft may be used to estimate when the Cape Cross NC has reached the distal end of the guiding catheter (depending on whether the approach is brachial or femoral).

Track the balloon catheter over the wire to cross the lesion, using the radiopaque markers to locate the balloon across the lesion.

Confirm that the balloon is correctly positioned.

Advance the dilatation catheter over the guide wire and into the stenosis or stent. Inflate the balloon (do not exceed 10 total inflations) to dilate the lesion using standard PTCA techniques. After each subsequent inflation, the distal blood flow should be assessed. If a significant restenosis persists, successive inflation may be required to resolve the stenosis. Maintain negative pressure on the balloon between inflations. Confirm the results with fluoroscopy.

Caution: Always advance and retract the Cape Cross NC fully deflated and always on the guidewire.

Notes:

- It is recommended that the guidewire and/or the balloon catheter remain across the lesion until the procedure is complete.
- Contrast media have different viscosities and may affect the inflation/deflation time. Use only diluted contrast medium to inflate the balloon (50% solution in

sterile heparinised saline).

- If two Cape Cross catheters are used simultaneously in a 6F guiding catheter (kissing balloon technique), position one guidewire and its associated balloon catheter at the first target site prior to inserting a second guidewire and its associated balloon catheter. When removing the catheters, remove one catheter and its associated guidewire completely prior to removing the second catheter and its associated guidewire.

8. CATHETER EXCHANGE PROCEDURE

Deflate the balloon fully by pulling negative pressure with an inflation device.

Loosen the haemostasis valve.

Hold the guidewire and the haemostasis valve in one hand, while grasping the Cape Cross NC catheter shaft in the other hand.

Maintain the guidewire's position in the coronary artery by holding the wire stationary and begin pulling the Cape Cross NC dilatation catheter out of the guiding catheter.

Note: during this exchange, monitoring of the guidewire may be done under fluoroscopy.

Pull the Cape Cross NC catheter until the opening in the guidewire lumen is reached.

Slide the remainder of the Cape Cross NC catheter along the guidewire until the tip of the Cape Cross NC catheter exits the haemostasis valve. Tighten the haemostasis valve onto the guidewire and hypotube to hold it securely in place.

Completely remove the Cape Cross NC catheter from the guidewire. Be careful not to twist, or rotate the Cape Cross NC balloon catheter around the guidewire.

Prepare and insert the next catheter (according to the manufacturer's Instructions for Use).

After the procedure, remove the guiding catheter from the vasculature.

9. ADVERSE EFFECTS

Possible adverse effects include, but are not limited to, the following:

- death,
- acute myocardial infarction,
- acute vessel closure,
- cardiac tamponade/pericardial effusion,
- coronary aneurysm or pseudoaneurysm,
- total occlusion of the coronary artery or bypass graft,
- coronary vessel dissection, perforation, rupture or injury,
- restenosis of the dilated vessel,
- hemorrhage, hematoma or other bleeding complications,
- unstable angina,
- arrhythmias, including ventricular fibrillation,
- drug reactions, allergic reaction to contrast medium,
- hypotension/hypertension,

- infection,
- coronary artery spasm,
- arteriovenous fistula,
- stroke, air embolism and embolization or fragmentation of thrombotic or atherosclerotic material,
- pain/anxiety,
- side branch occlusion,
- ischemia.

These complications may result in the need for percutaneous reintervention, coronary artery bypass graft surgery, or surgical repair or intervention.

10. STORAGE REQUIREMENTS

Use before the expiry date clearly indicated on the label.

The product is stored between 5°C and 25°C, away from light and moisture, outside the cleanroom. Long term exposure to high temperatures may reduce the shelf life. The product must not be exposed to extreme temperatures below 5°C and above 40°C.

11. WARRANTY AND LIMITATION OF LIABILITY

DISA Vascular guarantees that each of its products is treated with utmost care in development, selection of materials and components, manufacture, sterilisation and final testing before it is released for packaging and distribution. Through inexpert storing, handling and other manipulation, the catheter system can be damaged and its operability can be restricted. DISA Vascular accepts no liability or warranty for malfunction, breakdown and potential medical complications for patient and hospital staff, or for resulting damages arising from inexpert handling, operation, storing, from acts of God or from other external influences or manipulation.

DISA Vascular will replace any device, which, in its opinion, was defective at the time of shipment and if defects that were caused during manufacturing or packaging are immediately brought to the attention of DISA Vascular or its distributors. This warranty is exclusive and in lieu of all other warranties, whether expressed or implied, written or oral, including, but not limited to, any implied warranties of merchantability of fitness.

As a result of biological differences in individuals, no product is 100% effective under all circumstances. DISA Vascular has no control over the operation, inspection, maintenance, or use of its products after sale, and has no control over the selection of patients. Therefore, DISA Vascular and its distributors do not guarantee either a good effect or against a poor result following the use of any DISA Vascular product. DISA Vascular and its distributors shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of the device or from re-sterilisation or re-use of the product.

In the case of product complaint, an appropriate complaint protocol is to be obtained from DISA Vascular and returned, together with the product, to DISA Vascular.

12. COMPLIANCE

The compliance chart is based on data from in vitro testing at 37°C, rounded to reflect the nominal diameter at nominal pressure. All data is provided to 2 decimal places. Nominal pressure is 14 bar.

Rated burst pressure (RBP) is 22 bar. Do not inflate beyond RBP (rated burst pressure).

13. REFERENCES

The physician should consult recent literature on current medical practice on balloon dilatation, such as that published by the American College of Cardiology and the American Heart Association.

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The Cape Cross NC Catheter is a product of DISA Vascular 2015 (Pty) Ltd.