


CAPECROSS

R_x PTCA Balloon Catheter

Improved
product

BALLOON CROSSING PROFILE


Tip Entry Profile – 0.40mm (0.016") for all diameters



Balloon dimensions were the same for all products tested (3.0mm)

*Data on file at DISA Medinotec

TRACKABILITY



TECHNICAL INFORMATION

- Balloon Material – Polyamide
- Shaft Diameters (3.0 x 20mm) – 2.5F distal, 1.9F proximal
- Catheter Length – 142cm
- Minimum Guiding Catheter – 5F
- Nominal Pressure – 6Atm
- Rated Burst Pressure – 16Atm
- Kissing Balloon Compatible – 6F, up to 3.5mm

SIZES AND ORDER CODES

Balloon Diameters (mm)	Length (mm)					
	10	12	15	20	25	30
1.5	CCD1-1510	CCD1-1512	CCD1-1515	CCD1-1520	CCD1-1525	-
2.0	CCD1-2010	CCD1-2012	CCD1-2015	CCD1-2020	CCD1-2025	CCD1-2030
2.5	CCD1-2510	CCD1-2512	CCD1-2515	CCD1-2520	CCD1-2525	CCD1-2530
3.0	CCD1-3010	CCD1-3012	CCD1-3015	CCD1-3020	CCD1-3025	CCD1-3030
3.5	CCD1-3510	CCD1-3512	CCD1-3515	CCD1-3520	CCD1-3525	CCD1-3530
4.0	CCD1-4010	CCD1-4012	CCD1-4015	CCD1-4020	CCD1-4025	CCD1-4030



DISA Medinotec (Pty) Ltd
 Stand 171, Northlands Business Park, Bush Telegraph Street,
 Northriding, Johannesburg, 2169, South Africa
 T: +27 87 330 2301 www.disamedinotec.com



Please see the instructions for use for a complete listing of the indications, contra-indications, warnings and precautions.
 CapeCross PTCA Balloon Catheter
 A4 Brochure Issue 05

